
Info Sheet

High Pressure- Back Pressure Frame

Max Pressure	Temperature Range	Working Volume
1800 psig	30 - 330C	0.5 – 100 ml

Features:

- Designed to investigate any reaction or process that has homogeneous liquid feeds
- Equipped with a 6-port valve for on-line sampling
- A typical setup consists of Isco pumps, heated coiled tube reactor, and Back Pressure Cart
- Enables use of common solvents at elevated temperature and pressure (environmental and cost advantages)
- Cannot handle solids (ppt)